



10	MOG-65-15	Goniometer	1			0	
9	MOG-65-10	Goniometer	1		571.79	571.79	
8	20MP	Mounting Plate	1	2A12	133.15	133.15	
7	MOR-100-30	Rotary Stage	1		1553.262	1553.262	
6	15MP	Mounting Plate	1	2A12	278.96	278.96	
5	8MP	Mounting Plate	1	2A12	165.88	165.88	
4	GB/T70.1-2000	Screw M4x10	4	Q235	2.05	8.2	
3	15AB	Angle Bracket	1			0	
2	GB/T70.1-2000	Screw M6x10	31	Q235	5.30	164.3	
1	MOX-02-30	Linear Stage	3		1093.98	3281.94	
IndexID	PCPNo	PCPName	Amount	MaterialName	Weight	TWeight	Remark

标记	处数	更改文件号	签字	日期
设计	李国胜	主任设计		
校对		标准化	李国胜	
审核	朱秀雪	批准	张温暖	
工艺	李国胜	日期	2020/3/19	

6-Axis Positioning Stage

MOA-30

图样标记	件数	比例
S		1:5

共 1 张 张 数

Optics Focus